

1. install git
2. install and setup azcopy
3. install uv
4. open Git Bash and run

```
cd $HOME
```

```
git clone https://github.com/SurgeonDecisionIntelligence/quality-labeler.git
```

```
azcopy cp --recursive https://tissueconnectstorage.blob.core.windows.net/joshfisher/Irina/QC
```

```
cd QCLabeler
```

```
uv run download_qc.py lineage_092126.csv $HOME/Pictures/QC
```

1. After all the images download, open Git Bash and run

```
cd $HOME/quality-labeler
```

```
uv run quality-labeler label $HOME/Pictures/QC
```